

Non-viral targeted knock-in of a KRAS G12D specific TCR, CD8α/β, and chimeric cytokine receptor in the TRAC locus outperforms lentiviral-based engineering of T cells

¹Santosh Narayan, ¹Allison P Drain, ¹Nathaniel Swanson, ¹Nicholas Rouillard, ¹Tyler Warner, ¹Nicole Danek, ¹Ken Gareau, ¹Cheryl Black, ¹James Parsons, ¹Anthony Thomas, ¹Jinsheng Liang, ¹Luhua Shen, ¹Tanya Tetrault, ¹Vince Nguyen, ¹Iqraa Priyata, ¹Sarah Vidyasagar, ¹Joshua Francis, ¹Xingyue He, ²Patrick J. Browne, ²Rebecca C Lamothe, ²Meghan D Storlie, ²Gregory J Cost, ³Thomas M Schmitt, ^{3,4,5}Philip D Greenberg, ⁶Smita S Chandran, ⁶Christopher A Klebanoff, ¹Ankit Gupta, ¹Damien Hallet, ¹Gary Shapiro, ¹Kim Nguyen, and ¹Loic Vincent
¹Affiniti Therapeutics, Seattle, WA and Watertown, MA, USA. ²Metagenomi Inc., Emeryville, CA, 94608 USA., ³Fred Hutchinson Cancer Research Center, Seattle, WA, USA. ⁴University of Washington School of Medicine, Seattle, WA, USA. ⁵Departments of Immunology and Medicine, University of Washington, Seattle, WA, USA. ⁶Center for Cell Engineering, Memorial Sloan Kettering Cancer Center, New York, NY, USA.

Abstract

T cells engineered with T cell receptors (TCRs) enable targeting of clonally-expressed oncogenic driver mutations and have potential to induce durable responses in patients with solid tumors. Viral vectors, including lentiviral vectors (LVV), have been a standard modality to deliver transgenes for T cell therapies but are severely limiting in their manufacturing time, cost, and cargo size. In contrast, non-viral targeted gene knock-in (KI) overcomes these limitations, substantially reducing manufacturing complexity. Here, we compare primary human T cells engineered using LVV or KI processes to express a TCR recognizing a KRAS G12D mutant peptide presented on HLA-A*11:01, the CD8α/β coreceptor, and a chimeric interleukin receptor (ILR). We developed and optimized a non-viral manufacturing process that uses a novel CRISPR-Cas12a system to knock-in the transgene cassette within the *TRAC* locus and to simultaneously knock out the endogenous TCR¹. Our non-viral KI platform edits primary T cells with high efficiency, and we show that KI engineered TCR T cells performed equivalently or better in functional assays relative to LVV-engineered cells. Together, these data support the utility of a non-viral gene KI approach and its planned incorporation into clinical development.

Overview of Transgenes

Autologous primary human CD4+ and CD8+ T cells were engineered to express:

- A high avidity HLA-A*11:01-restricted TCR specific for the *KRAS G12D* mutation,
- CD8 α/β co-receptor that drives a coordinated CD8/CD4 T cell response by allowing for CD4 stimulation that promotes CD8+ T cell functional persistence,
- Interleukin Receptor, ILR, a fusion protein that promotes anti-tumor activity through increased T cell proliferation and survival.




Fig 1: Lentiviral delivery is limited by cargo size




Fig 2: Non-viral Knock-in can achieve high transgene integration frequency even with large transgenes




Fig 3: KI TCR T cells show improved tetramer binding




Fig 4: KI TCR T cells show high functional avidity to the KRAS G12D peptide




Fig 5: KI TCR T cells show robust cytotoxicity in vitro




Fig 6: KI cell show superior anti-tumor activity in a xenograft model in vivo




Fig 7: EF-1α and Promoter-less construct designs for KI




Fig 8: EF-1α promoter drives stronger transgene expression




Fig 9: EF-1α promoter resulted in higher functional avidity





Fig 10: EF-1α promoter drives superior anti-tumor activity in an in vivo model



Summary

- > Non-viral targeted gene knock-in achieved high KI efficiency and cell expansion in primary CD4/CD8 T cells.
- > KI engineered cells show higher transgene expression and superior *in vitro* and *in vivo* activity compared to cells engineered with LVV despite a lower transgene copy number per cell
- > EF-1α promoter drives higher transgene expression than the endogenous *TRAC* promoter leading to increased avidity and functionality of KI cells
- > The KI process has been established for the AFNT-212 product poised to enter clinical testing in 2024.

References

1. Goltsman, D.S.A. et al. Novel Type V-A CRISPR Effectors Are Active Nucleases with Expanded Targeting Capabilities, *CRISPR J.*, 2020.
2. Schober, K. et al. Orthotopic replacement of T-cell receptor α- and β-chains with preservation of near-physiological T-cell function, *Nature Biomed. Eng.*, 2019.

See SITC Abstract #355 for AFNT-212 product overview